



TEOSYAL® RANGE - INDICATIONS AND SUMMARY TABLE

	TEOSYAL® Meso	TEOSYAL® First Lines	TEOSYAL® Global Action	TEOSYAL® Touch-Up	TEOSYAL® Deep Lines	TEOSYAL® Kiss	TEOSYAL® Ultra Deep	TEOSYAL® Ultimate
Concentration	15 mg/g	20 mg/g	25 mg/g	25 mg/g	25 mg/g	25 mg/g	25 mg/g	22 mg/g
Volume-creating capability	○○○○○	●○○○○	●●○○○	●●●○○	●●●○○	●●●○○	●●●○○	●●●●●
Needle	30G 1/2	30G 1/2	30G 1/2	30G 1/2	27G 1/2	27G 1/2	26G 5/8	21G 1 1/2
Indication	Rehydration of the dermis with treatment protocol.	Lines and superficial crow's feet and peribuccal wrinkles.	All moderately pronounced wrinkles. Glabella wrinkles, moderately pronounced nasogenian folds and peribuccal wrinkles.	All moderately pronounced wrinkles, ideal for touch-ups.	Highly pronounced wrinkles. Nasogenian folds and deep vertical lip lines. Moderately pronounced lower eyelid depression.	Lip contour and volume.	Creates and restores volumes, facial contours and lower eyelids. Very deep wrinkles in thick skin.	Cheeks, upper cheeks, chin, facial contours.
Skin lesions								
Blister pack	2 x 1ml	2 x 0.7ml	2 x 1ml	2 x 0.5ml	2 x 1ml	2 x 1ml	2 x 1ml	1 x 3ml
Estimated duration before top-up*	2 - 3 months	6 months on average	6 - 9 months on average	6 - 9 months on average	6 - 9 months on average	6 - 9 months on average	9 - 12 months on average	12-18 months on average

*The average lifespan depends on several factors: patient's skin type, severity of the wrinkle to be corrected, type of injection (superficial, medium or deep dermal) and the volume injected.

LABORATOIRES
TEOXANE
GENEVA

Les Charmilles - Rue de Lyon, 105 - CH 1203 Genève
Tél. : +41 (0) 22 344 96 36 - Fax : +41 (0) 22 340 29 33
info@teoxane.com

THE BEST OF HYALURONIC ACID



TEOSYAL Ultimate
Intense
volume

LABORATOIRES
TEOXANE
GENEVA